

DRE Milano E20

Power Procedure Chair

Advanced features and elegant engineering

The DRE Milano E20 is an excellent alternative to traditional procedure chairs. It offers innovative features and modern styling at a competitive price. The multi-specialty general procedure chair is available with options suited for most office-based medical procedures.

The Milano E20 is comfortable for both the operator and patient: it is wider than traditional procedure chairs and offers easy access to people with mobility difficulties, thanks to the chair's adjustable height (lowers to 20.5"/52 cm).

Designed for comfort, the Milano E20 can be used for lengthy procedures. Plus, its physician-configured adjustments provide excellent visibility within the operative field.

Features

- ✓ 180° chair rotation: Swivel feature provides for efficient access to multiple patient areas.
- ✓ Tapered head allows for easy access to facial areas for the practitioner.
- ✓ Foot pedal adjusts the height, back and level positioning.
- ✓ Lumbar support mattress: Ergonomic mattress design allows for lumbar-back support, patient head comfort and adjustable foot support.
- ✓ Convenient "auto return" function restores the chair to original position and saves operator time. This allows the patient to easily get on and off the chair.
- ✓ Footrest is easily extended via a manual control knob.
- ✓ Sleek system with a small footprint.
- ✓ Stitchless mattress/covering allows for easy cleanup and sterilization.



374-01-2017-01

SPECIFICATIONS



Chair Weight: 166 lbs. (75kg)

Maximum Capacity Load:
440 lbs. (200kg)



Voltage: 100-110 - 230V


Frequency: 50/60 Hz

Maximum Consumption:
100W



Phases: L+N

Color Options



Optional Accessories

- ✓ Standard U.S armboard rails for a variety of accessories
- ✓ Roll holder
- ✓ Lamp adaptor
- ✓ Heel support
- ✓ Metal tray with support
- ✓ Jointed thigh support
- ✓ For bariatric patients: Can be manufactured with 4 additional inches on the seat, taking it from 23" to 27".

