

Portable Halogen Headlights

DRE Xavier-C2

Maximize the view during examinations with the versatile Xavier-C2 Portable Halogen Headlight System.

- » Comfortable, adjustable Standard or Sport Headband.
- » Precision iris adjustment provides a variable spot size.
- » Lens and mirror system delivers clean, even light into the field.
- » Brilliant 11.5 cm spot at 40 cm working distance.
- » Autoclavable positioning bar provides sterile adjustment.
- » Pivot points allow for exceptional coaxial alignment.
- » Component design for low cost and simple maintenance.
- » Brilliant 10-watt halogen lamp provides 100 hours of use.
- » In-line on/off switch provides easy access control.



DRE Xavier-MH

The Xavier-MH is the solution when you need to deliver intense, shadow-free light into any body cavity. Designed for otolaryngology, it provides comfort, intensity and versatility.

- » Comfortable, adjustable headband.
- » Swivel technology allows the optic to be precisely positioned over the dominant eye for true monocular viewing.
- » Brilliant 11.5 cm spot at 40 cm working distance.
- » Three Mag-Fluoride coated lenses deliver focus and intensity.
- » Autoclavable positioning bar provides sterile adjustability.
- » Component design for low cost, simple maintenance.
- » Brilliant 10 watt Halogen Lamp provides 100 hours of use.
- » In-line on/off switch provides easy access control.



Battery-operated – no light source required!

Two-hour, rechargeable belt-clip battery for mobile power.

110v - 240v transformer delivers uninterrupted use.

110v or optional 240v charger included with each system.